

# Gasperini "Shiny Tuscany" thermographic powders

# **Neutral Gloss**

Gloss 7 - extra coarse

Gloss 9 - coarse

Gloss 11 - medium coarse

Gloss 14 - medium fine

Gloss 18 - extra fine

# **Neutral Matt**

Matt 14 - medium fine

Matt 18 - extra fine

#### **Pearlescent**

Pearl Satin Gold 14

Pearl Satin Copper 14

Pearl Satin Brass 14

Pearl Satin Silver 14

Pearl Satin Blue 14

Pearl Satin Green 14

Pearl Satin Red 14

Pearl Satin Violet 14

Pearl Satin Orange 14

#### **Fluorescent**

Fluorescent Pink 14

Fluorescent Green 14

Fluorescent Orange 14

web:

e-mail:

www.gasperini.it

info@gasperini.it

C.F.: 03682920487

P.IVA/VAT No.: IT00317750974

Fluorescent Blue 14

# **Opaque**

Opaque Titanium White 14 Opaque Porcelain White 14

# **Special**

White Puff-Up 11

Glow in the Dark

# **Granulation sizes**

granulation	use
07	Extra coarse. Used for large and solid image areas.
09	Coarse. Generally used for large image areas and types between 36 and 72 points.
11	Medium-coarse. Generally used for large/medium image areas and types between 27 and 42 points.
14	Medium-fine. It is the most used granulation for its excellent versatility. Suitable for types between 11 and 30 points.
18	Extra fine. It is the classic granulation for fine and italic types between 6 and 14 points.